## Fingertip Pulse Oximeter Instructions for Use SH0-1002

### 1. Safety information

- Do not use the pulse oximeter in an MRI or CT environment.
  - Do not use the pulse oximeter in situations where alarms are required. The device has no alarms.
  - Explosion hazard: Do not use the pulse oximeter in an explosive atmosphere.
  - The pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
  - Check the pulse oximeter sensor application site frequently to determine the positioning of the sensor and circulation and skin sensitivity of the patient.
  - Do not stretch the adhesive tape while applying the pulse oximeter sensor. This may cause inaccurate readings or skin blisters.
  - Before use, read the manual carefully.
  - The pulse oximeter has no SpO2 alarms; it is not for continuous monitoring.
  - Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status, and correct alignment at least every 4 hours.
  - Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.
  - Significant levels of dysfunctional hemoglobins (such as carboxy-hemoglobin or methemoglobin) may affect the readings.
  - Intravascular dyes such as indocyanine green or methylene blue
  - SpO2 measurements may be adversely affected in the presence of high ambient light or direct sunlight. Shield the sensor area (with a surgical towel, for example) if necessary.
  - Excessive patient movement may cause inaccurate readings.
  - High-frequency electrosurgical interference may cause inaccurate readings.
  - Venous pulsations may cause inaccurate readings.
  - Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
  - If the patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia, this may cause inaccurate readings.
  - Do not use if the patient is in cardiac arrest or is in shock.
  - Nail varnish or false fingernails may cause inaccurate SpO2 readings.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.


### 2. Introduction

The principle of the Oximeter is as follows: A mathematical formula is established making use of the Lambert Beer Law on the Spectrum Absorption Characteristics of Reductive hemoglobin (RHb) and Oxyhemoglobin (HbO2) in glow and near-infrared zones. Operation principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelength of lights (660nm glow and 905 nm near infrared light) can be focused onto a human nail tip through a finger clamping sensor. A measured signal

obtained by a photosensitive element will be shown on the Oximeter's display through a process in electronic circuits and a microprocessor.

### Diagram of Operation Principle

- Red and Infrared-ray Emission Tube
- Red and Infrared-ray Receipt Tube


The Fingertip PULSE OXIMETER is a portable non-invasive, spot-checkfor oxygen saturation of arterial hemoglobin (SpO2) and pulse rate of adult and pediatric patients at home, and in hospital (including clinical use in internist/surgery, anesthesia, intensive care etc). It is not for continuous monitoring.

### 3. Features

- Operation of the product is a simple and convenient, easy operation without pressing a button.
- The product is small in volume, light in weight and convenient to carry.
- Power consumption of the product is low and the two AAA batteries can be operated continuously for 30 hours.
- A low voltage warning will be indicated in a window when battery voltage is so low that normal operation of the oximeter might be influenced.
- The product will automatically be powered off when no signal is in the product for longer than 5 seconds.
- No maintenance is needed except to change batteries.

### 4. Product introduction

### 5. Operation

- Install two AAA batteries into the battery cassette before closing its cover.
- Auto power on as follows.

- Insert one finger fully into the rubber hole of the Oximeter.
- Finger and body should not tremble during measurement.
- Read corresponding data on the display screen directly.
- If the finger leaves the rubber hole for 5 seconds, it can auto off.

**Note: Please use medical alcohol to clean the rubber inside the Oximeter. Also, clean the test finger with alcohol before and after each test. (The rubber inside of the Oximeter is medical rubber, which has no toxins, and is not harmful to the skin).**

### 6. Specifications

#### 6.1 Performance

- Display Type: OLED
- SpO2:
  - Measurement range: 35%-100%
  - Accuracy: 70%-99%, ±2% ; ≤69%, no definition.
- Pulse Rate:
  - Measurement range: 30-250BPM
  - Accuracy: 30-250BPM, ±%
  - Pulse rate display model: Bargraph Indicator

#### 6.2 Electrical specifications

- Power Requirements:
- Two AAA alkaline Batteries
  - Power consumption: Less than 30mA
  - Low power indication: The low voltage lamp lights up when the voltage less than 2.6V (error±0.2V)
  - Battery Life: Two AAA 1.5V, 600mAh alkaline batteries can be continuously operated as long as 30 hours.

#### 6.3 Environmental conditions

- .Environment Requirements:.
- Operation Temperature: 5-40°C Storage Temperature: -20-70°C
  - Ambient Humidity: Work: 15%-90% Storage: 15%-93

#### 6.4 Physical characteristics

- .Dimensions:
- L×W×H : 57mm×34mm×30.9mm(SHO1001)
  - 57mm×36mm×32mm(SHO2001)
  - 57mm×31.5mm×30.5mm(SHO3001N)
  - 57mm×31mm×31mm(SHO4001)

#### 6.5 Standards

- EMC of this product complies with IEC60601-1-2 standard.
- Measurement Performance in Low Perfusion Condition: This pulse oximeter gives normal results whenever the wavelength of the pulse simulation (FLUKE-INDEX 2XL) is 0.4% 3. Interference Resistance
- Capacity against Ambient Light: Device works normally when light interference is produced by FLUKE-INDEX 2XL.

#### 6.6 Markings

Symbol	Definition
	Type B applied part.
	Attention, refer to the manual before using.
%SpO2	Oxygen saturation
PRbpm	Heart rate (BPM)
	Low power indication
	No SpO2 Alarm
SN	Serial No
	Do not discard in garbage bin

### 7. Maintenance and Storage

- Replace the batteries in time when low voltage lamp is lighted.
- Clean surface of the fingertip oximeter before it is used in diagnosis of other patients.
- Remove the batteries inside the battery case if the Oximeter will not be operated for a long time.
- It is best to keep the product in a place where ambient temperatures is -10°C ~ 40°C (14 ~ 104°F) and relative humidity is between 15% ~ 80%.
- It is recommended that the product should be kept in a dry environment at all times. A wet environment might affect its lifetime and even might damage the product.
- Please follow the law of the local government to deal with used batteries.

### 8. Trouble-shooting and maintenance

Problems	Possible Reason	Solution
SpO2 or PR can not be shown normally	Finger is not inserted correctly	Retry by inserting the finger
SpO2 or PR is shown	1.Finger might not be inserted deeply enough. 2. Finger is trembling or patient's body is in movement status.	1.Retry by inserting the finger 2. Try not to move
The Oximeter can not be powered on	1. Power of batteries might be inadequate or not be there at all 2.Batteries might be installed incorrectly 3.The Oximeter might be damaged	1.Please replace batteries 2.Please reinstall the batteries 3.Please contact local customer service centre
Indication lamps are suddenly off	1.The product is automatically powered off when no signal is detected longer than 5 seconds 2. The batteries are running down.	1.This is normal 2.Replace the batteries

### 9. Warranty

- The company guarantees the pulse oximeter at the time of its original purchase and for the subsequent period of two year. Spare parts guarantee is six months. The warranty does not cover the following:
- If the device series number label is torn off or cannot be read.
  - Damage to the device resulting from misconnection with other devices.
  - Damage to the device resulting from accidents.
  - Changes performed by users without the prior written authorization of the company.

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Technonology Ltd.