Invacare® Pausa (L865)  
 Invacare® Kauma (L872)  
 Invacare® Areste (L110)  
 Invacare® Satio (L871)

|                      |    |
|----------------------|----|
| Guide utilisateur    | FR |
| User manual          | EN |
| Bedienungsanleitung  | DE |
| Manual del usuario   | ES |
| Manuale d'uso        | IT |
| Gebruiksaanwijzing   | NL |
| Manual de Utilização | PT |
| Bruksanvisning       | SV |
| Brugermanual         | DK |
| Käyttäjän käsikirja  | FI |
| Bruker manual        | NO |

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## 2. ENGLISH

You have just purchased an INVACARE® product. Our products are designed and manufactured to the most rigorous quality and safety standards and conform to the European directives (MDD 93/42/EEC). They are guaranteed for one year. Certain products may have a longer guarantee; please check with your dealer.

**The following operating and maintenance instructions must be observed:**

- Check that the product has been delivered in perfect condition before you use it.
- If you have any doubts about the restrictions on use, the adjustment or maintenance of the product, please contact your dealer.
- Do not exceed the maximum operating load  or the maximum load  shown on the product label.
- Always ensure that the product is suitable for your requirements; in case of doubt, check with a health professional.
- Worn, damaged or missing components must be replaced immediately by an approved repairer.
- Do not leave or store the product out of doors.
- Periodically check the stability and rigidity of the product and, where applicable, that the screws are correctly tightened and the rubber end caps in good condition.

### L865 BED TABLE

- Fix the 2 feet on castors under the base of the "L" using 2 hex screws and insert the adjusting screw into the upper section of the "L".
- Slide the square part of the tray support into the vertical section of the "L"; the tray pivot should sit above the castor base.
- Unscrew the tray pivot adjusting screw and slide the pivot into position below the tray.
- Adjust the tray to the desired inclination and tighten the tray pivot adjusting screw; set the height of the tray and tighten the main adjusting screw.

### L872 BED TABLE

Warning: the height adjustment system is controlled by a spring. Follow the instructions in this paragraph carefully, in the order indicated, when fitting the table.

#### Fixing the upright on the castor base

- Insert the two large screws in the two-hole plate, and then insert the two small cylinders. Place the whole assembly on the castor base and screw it onto the upright; the tray pivot should sit above the castor base.

#### Fixing the tray

- Slide the tray onto its support and fix it using two 8x15 screws.

#### Fitting the adjustment system

- Insert the black lever, with the rounded section uppermost, into the horizontal section and lock it in place (using a 5x50mm bolt and nut).
- Carefully unscrew the small black screw above the red label and replace it with the small grey screw supplied with the table. Never remove the screw covered by the red label.
- Remove the protective film from the tray.

#### Use

- To raise the tray, pull it upwards.
- To lower the tray, press the black lever.

### L110 SHOULDER SUPPORT

Adjust the inclination and check that the shoulder support is properly locked.

### L871 BED CRADLE

Position the cradle above the patient. Arrange the covers over the cradle.

### MAINTENANCE

Remove any dirt or dust using a soft, damp cloth. Use soap and water to remove stubborn stains. Dry with a soft cloth. For disinfection in hospital environments, please follow the hygiene specialist's recommendations.

### USER'S MAXIMUM WEIGHT - MAXIMUM LOAD

Please see the table at the end of this document.

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